

LIFEPAK® CR2

defibrillator cellular


Features

- Easy to use¹
- Cellular connectivity
- Self-monitoring
- Location services¹
- Child mode
- Bilingual

You can't predict when an emergency strikes, but you can be prepared.

Sudden cardiac arrest (SCA) can happen to anyone—anywhere. Immediate treatment is vital. A victim's chance of survival dramatically decreases for every minute without treatment². That's why public access defibrillators are so important.

- **Metronome and CPR coaching**
Sets an effective pace and audibly guides users³.
- **QUIK-STEP® electrodes**
Peel directly off the base for faster placement¹.
- **Child mode**
Child mode delivers lower energy levels appropriate for young children without having to change electrodes.
- **ClearVoice™ technology**
Detects background noise and adjust tones and voice prompts to ensure they can be heard clearly in noisy environments.
- **Cellular connectivity**
Device will use a 4G cellular connection and/or Wi-Fi to communicate with LIFELINKcentral™.
- **LIFELINKcentral AED Program manager**
Monitor AED programs by tracking AED status, sending patient data to hospitals, and giving device specific location information.



References:

- 1 Physio-Control Internal Semi-Automatic AED Comparison Usability Study, August 2016.
- 2 Graham R, McCoy M, Schultz A. Strategies to Improve Cardiac Arrest Survival, A Time to Act. Institute of Medicine Report, 2015.
- 3 Mayo Foundation for Medical Education and Research. (2024, April 30). Cardiopulmonary resuscitation (CPR): First aid. Mayo Clinic. <https://www.mayoclinic.org/first-aid/first-aid-cpr/basics/art-20056600#:~:text=Use%20your%20entire%20body%20weight,spring%20back%20after%20each%20push.>

This document is intended solely for the use of healthcare professionals. A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate Stryker's product offerings. A healthcare professional must always refer to operating instructions for complete directions for use indications, contraindications, warnings, cautions, and potential adverse events, before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Specifications subject to change without notice.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: ClearVoice, cprINSIGHT, LIFEPAK, LIFEPAK TOUGH, LIFELINKcentral, QUIK-STEP, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.



Manufactured by:
Physio-Control, Inc.
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
stryker.com/emergencycare

Distributed in Canada by:
Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323